

PROCESS OR MATERIALS REVIEW REQUEST FORM

Request Title: Centrifugation of polymer-metal complex

Requestor: Dr. Ashutosh Divekar

Date: 06/27/2023

Non-SNF staff, please complete the following:

Badger ID: ashutosh.divekar@evoloh.com

Company name/PI: Evoloh inc.

Send completed form to: snf-promcommittee@lists.stanford.edu

1. The chemical or material.

The chemical/material of interest is a solution with insoluble solids. It is an aqueous solution containing 33.3% dimethyl sulfoxide (DMSO), 65% H₂O, and 1% polymer-metal complex.

The solution doesn't have an MSDS as it will be a mixture post-synthesis. However, for synthesis, the following chemicals will be used: dimethyl sulfoxide (solvent), hydrocarbon polymer (polymer), porphorin (cross-linker), cobalt acetate (metal salt).

The solvent DMSO is flammable; however, flammability is insignificant at 33% by volume. Hydrocarbon polymer-metal complex is not hazardous and is just ~1 % of the total solution.

I only need the centrifuge capabilities, so I will not store the chemical post-centrifugation. But, the storage group would be "G" as it is mostly an aqueous solution.

2. Vendor/manufacturer info:

Solutions will be generated in-house at Evoloh inc. 3503 Edison Way, Menlo Park, CA 94025.
720-717-8735.

3. Reason for request:

Evoloh does not have a centrifuge in-house, so it needs SNF capabilities.

4. Process Flow:

The solution will be stored into compatible (polypropylene) centrifuge vials. The vials will be sealed with parafilm and placed into another secondary container. The secondary container will be brought into the Exfab Mavericks 155 lab in a secondary container. The vials will then be placed next to the centrifuge located in one of the hoods. The vials will then be centrifuged for a set duration of time. Post-centrifugation, the vials will be placed inside the

PROCESS OR MATERIALS REVIEW REQUEST FORM

secondary container and directly taken out of the facility.

5. **Amount and form.**

The solution will be brought in solids suspended in the liquid. The solution will be in centrifuge vials tightly sealed with parafilm. The solution will contain ~1% solids wt/vol. The total volume will be 300 ml which will be split into multiple centrifuge vials as per the size compatibility with the centrifuge equipment (for e.g., 6 vials of 50 ml each). The solution will only be inside the Ex-fab Mavericks 155.

6. **Storage:**

The solution will not be stored at SNF.

7. **Disposal:** The solution will be disposed of internally at Evoloh and won't need to be disposed of at Stanford.

To be completed by the PROM COMMITTEE:

Special handling requirements:

On-going approval or staff interaction requirements:

Approved By/Date: